

# Lovibond<sup>®</sup> Colour Measurement Tintometer<sup>®</sup> Group



## Liquid Standard Platinum-Cobalt Colour 15 - 60 ml

Part Number: 135049



Reference liquid standards are assigned a nominal value and Lovibond Tintometer always endeavours to match this value as closely as possible. On occasion, it may not be possible to match the requested value exactly. The value achieved and the expected performance tolerances will, of course, always be reported.

### Liquid Standard Platinum-Cobalt Colour 15 - 60 ml

60 ml Bottle certified under ISO 9001 quality system with a 12-month shelf life.

Full traceability to internationally recognised standards

Each bottle supplied with full certification including MSDS