

Lovibond® Colour Measurement Tintometer® Group



Hazen Colour Filter Set - EC 2000

Part Number: 135159



Conformance filters are assigned a nominal value and the company always endeavours to match this value as closely as possible. On occasion, because of the unique manner of the glass making process, it may not be possible to match the requested values exactly. The values achieved and the performance tolerances expected will, of course, always be reported.

Hazen Colour Filter Set - EC 2000

Glass conformance filters for checking Pt-Co scale on EC 2000 PT-Co.

Technical Data

Glass Standard Values	Approximate values
	- 5
	- 20
	- 50
	- 100
	- 300

Filters 5

Delivery Scope

- 5, 20, 50, 100, 300 filters
- Certificate of Conformity