

Lovibond® Water Testing

Tintometer® Group



Legionella Field Test Kit

The World's Fastest *Legionella* Test



Legionella is the causative agent of Legionnaire's disease, a potentially fatal form of pneumonia. Guidelines exist across the world for the control of *Legionella*.

In 2002, the Health and Safety Executive released guidelines for the control of *Legionella*, known as the L8 guidelines. These guidelines detail the key parameters that should be tested on a weekly and monthly basis to minimise the risk of *Legionella* proliferating to harmful levels in open water system.

The Lovibond® *Legionella* "L8" Test Kits contain the recommended tests for both weekly and monthly monitoring.

- Fast, easy and accurate
- Results in as little as 25 minutes!
- Confirms performance of *Legionella* prevention and treatment protocols
- Test Kits available for detecting both waterborne and biofilm *Legionella*

Part Number: 56B006001

Legionella Field Test Kit

The *Legionella* Field Test Kit is the most basic of the *Legionella* Detection test kits. The test is easy and fast — no reagents or interpretation of results required!

With results in 25 minutes, the test allows for immediate understanding and response to the water conditions. Regular use will increase confidence in your *Legionella* control procedures and reduce risk of outbreaks and litigation!

Users simply collect a water sample in the provided collection bag, and using the provided exact volume disposable pipet, dispense 100 µL of sample onto the test strip. After the required reaction time, the test line will turn red if results are positive.

Limit of Detection
100,000 CFU/l

No. of Tests
10

Delivery Scope

- 10 Test Strips with Dispensing Pipettes
- 10 Collection bags
- Instructions Manual
- Carrying Case