

Lovibond® Water Testing

Tintometer® Group



Water Safety Kit Microbiology Duo

Mobile microbiological test equipment with 2 incubators



Disaster situation or potable water monitoring in demanding situations: You need certainty and reliability for test results in everywhere. The Water Safety Kit Microbiology Duo contains all the equipment necessary to perform the microbiological analysis component for safe drinking water. Due to the heating and cooling function of the DI 20 incubator you can test different sorts of bacteria - with both incubators simultaneously! Perform filtration with sterilisable equipment including sampler, filter flask & vacuum pump on-site - sterilisation equipment not included. With the optional backpack you can take the kit everywhere with you.

- Testing coliforms and E. coli bacteria according to WHO standards
- Testing different forms of bacteria simultaneously with 2 incubators
- DI 20 with heating and cooling as standard

Part Number: 56K681255

Industry

Industries Others | NGO

Application

Drinking Water Treatment

Water Safety Kit Microbiology Duo

In one hand or on the back: This test kit accompanies you in every situation with the complete microbiological test equipment. Test different sorts of bacteria directly on-site - with both incubators even simultaneously! Be sure of the quality of potable water everywhere, even when there is no laboratory nearby. The DI 20 incubator with its unique capability of heating and cooling as standard ensures reliable results.

User

Non Profit Organizations (NGO's) | Medical research and development | Institutes, Universities, Schools | World Health Organization (WHO)

Technical Data

Dimensions 487 x 229 x 386 mm

Weight with Packaging 18.75 kg

Delivery Scope

- Marker Pen
- 2.5" Bailer with Rope
- Filter Flask, 500ml
- Tweezers
- Pump, 35ml stroke
- 5 Universal Container and spoon, 30 ml
- DB 513 Shaker Tube with lid, 100 ml
- 150 EZ Pack Filters
- 5 Pipette Disposable, 1.5 ml (EDT222)
- 2 DI 20 Incubator
- 200 Media Pads and dispenser
- Incubator Battery Lead/Acid
- Filtration Device
- 50 Petri Dish for DI 20 Incubator
- Lead Set IEC C13, (3 pin)
- Membrane Lauryl Sulphate Broth, 40 g
- Instruction manual for Water Safety Kit Microbiology
- Instruction manual for DI 20 Incubator
- Case

Accessories

Title	Part Number
Power cable set	136300
Incubator Battery, Lead/Acid	192345
Bottles	363590
Universal container with spoon and lid, 30 ml	366239
Plastic shaker tube, graduated 100 ml, with lid	385130
Petri dishes, 10 pc.	400855
Backpack for Water Safety Kits	56A01400
Media Pads and Dispenser (Pack of 200)	56A016330
Bailer with rope, 2.5"	56A018501
EZ Filters 0.45 µm, 47 mm (pack of 150)	56A023801
Tweezer, 120 mm	56A604201
Filtration Device	56B000710
D004 TTC/MAC Dipslides (Macconkey) (pack of 10)	56B010410
D005 TTC/Chromogenic Dipslides (pack of 10)	56B010510
Membrane Lauryl Sulphate Broth, 38.1 gr.	56P069740
Pump (35 ml stroke)	56X610301
Pipette, disposable, 1.5 ml	BI-K10060
Marker Pen	BI-K14087

Tintometer GmbH

Lovibond® Water Testing
Schleefstraße 8-12
44287 Dortmund
Tel.: +49 (0)231/94510-0
sales@lovibond.com
www.lovibond.com
Germany

The Tintometer Limited

Lovibond House
Sun Rise Way
Amesbury, SP4 7GR
Tel.: +44 (0)1980 664800
Fax: +44 (0)1980 625412
sales@lovibond.uk
www.lovibond.com
UK

Tintometer China

9F, SOHO II C.
No.9 Guanghualu,
Chaoyang District,
Beijing, 100020
Customer Care China Tel.: 4009021628
Tel.: +86 10 85251111 Ext. 330
Fax: +86 10 85251001
chinaoffice@tintometer.com
www.lovibond.com
China

Tintometer South East Asia

Unit B-3-12, BBT One Boulevard,
Lebuh Nilam 2, Bandar Bukit Tinggi,
Klang, 41200, Selangor D.E
Tel.: +60 (0)3 3325 2285/6
Fax: +60 (0)3 3325 2287
lovibond.asia@tintometer.com
www.lovibond.com
Malaysia

Tintometer Brazil

Caixa Postal: 271
CEP: 13201-970
Jundiaí – SP
Tel.: +55 (11) 3230-6410
sales@lovibond.us
www.lovibond.com.br
Brazil

Tintometer Inc.

6456 Parkland Drive
Sarasota, FL 34243
Tel: 941.756.6410
Fax: 941.727.9654
sales@lovibond.us
www.lovibond.us
USA

Tintometer India Pvt. Ltd.

Door No: 7-2-C-14, 2nd, 3rd & 4th Floor
Sanathnagar Industrial Estate,
Hyderabad, 500018
Telangana
Tel: +91 (0) 40 23883300
Toll Free: 1 800 599 3891/ 3892
indiaoffice@lovibond.in
www.lovibondwater.in
India

Tintometer Spain

Postbox: 24047
08080 Barcelona
Tel.: +34 661 606 770
sales@tintometer.es
www.lovibond.com
Spain